DENTSPLY
DETREY

AquaCem[®]

Directions for Use _____	2
Gebrauchsanweisung _____	5
Mode d'emploi _____	8
Istruzioni per l'uso _____	11
Instrucciones de uso _____	14
Bruksanvisning _____	17

AquaCem®

AquaCem® is a light-yellow, translucent glass-ionomer luting material consisting of a blend of aluminosilicate glass and polyacrylic acid. The powder is mixed with distilled water to produce a luting material which adheres to dentine and enamel producing tightly sealed cementations.

AquaCem complies with the requirements of ISO 9917:1991 for polyalkenoate luting cements.

COMPOSITION

1000mg powder contain:	Calcium-sodium-fluoro-phosphoro-aluminium-silicate (21:2:15:2:14:46)	744.14 mg
	Polyacrylic acid	247.68 mg
	Tartaric acid	8.10 mg
	Yellow Ferric Oxide (E 172)	

INDICATIONS

Cementation of crowns, inlays, bridges, and orthodontic bands.

CONTRAINDICATIONS

- Direct or indirect pulp-capping.
- Cementation of extensive bridgework when the time needed for application of the cement and positioning of the bridge exceeds the working time (2 min 30 sec) of AquaCem.

SIDE EFFECTS

Pulpal irritation may occur in rare cases with glass-ionomer luting cements.

NEGATIVE INTERACTIONS WITH DENTAL MATERIALS

None reported.

DOSAGE

1 scoop powder : 2 drops water.¹

Fill water dispenser with distilled or de-ionized water. To dispense water, hold the bottle vertically above pad and squeeze gently.

Invert powder bottle before use to fluff the powder. Overfill powder scoop and strike off the excess with the bottle insert leaving a flat surface level with the edge of the scoop.

MIXING

Mix on the pad provided or on a glass slab.² Special spatulas are not required.

Divide the powder into two halves. Add one half to the water. As soon as the first half is wetted, incorporate the second, then mix thoroughly. The measures are intended as a guide; the mix should have a smooth creamy consistency.

It is essential that the powder is incorporated **very quickly** into the water to ensure the optimum mix ratio. For all normal uses the mixing time should not exceed 15 seconds.²

WORKING TIME²

The working time from start of mix is approximately 2 min 30 sec.

SETTING TIME³

The setting time from end of mix is 3 1/2 to 5 minutes.

APPLICATION

Dentine Pre-Treatment of the Preparation

Pulp Protection

Dentine in close proximity to the pulp should have been covered with a calcium hydroxide liner (Dycal®). The preparation is cleaned (e.g. prophylactic paste and/or H₂O₂ 3%), washed, and dried, but not dehydrated.

Application

Apply cement onto the fitting surface of the restoration/orthodontic band and/or the prepared tooth immediately upon completion of the mix whilst appearance is glossy. Position the restoration/orthodontic band and apply adequate pressure to seat it. Avoid contamination with water/saliva during application and setting. It should, however, be noted that the cement remains subject to dehydration if left isolated for longer than 10 minutes from start of mix. This means that water balance must be maintained by releasing the cement to the oral environment within this time range.

Finishing

Excess material may be removed at once or after completion of setting time (see above). Do not touch the cement when it is in the intermediate, rubbery condition. For removal of surplus material, appropriate hand instruments, like a dental probe, and dental floss are used.

CLEANING OF INSTRUMENTS

Glass-ionomers adhere to instruments and glass slabs. **Immediately after use**, either clean instruments with water or drop in disinfectant or water to facilitate subsequent cleaning.

SPECIAL NOTES

For dental use only.

Keep away from children.

Avoid contact with eyes. In case of contact, rinse eyes immediately with plenty of water and seek medical attention.

STORAGE AND SHELF-LIFE

Not to be stored at temperatures exceeding 25 °C.

Replace cap tightly after use.

Under these conditions, AquaCem has a minimum shelf life of three years.

Humidity can adversely affect the properties of glass-ionomer cements.

The desiccant crystals will protect the powder by absorbing the moisture in the bottle, but they are no longer effective after the colour has turned to pink

¹ The same powder:liquid ratio (3.3 g to 1 g) is proposed for testing purposes according to ISO 9917:1991.

² Mixing and application of the material should be done at normal room temperature. Higher temperatures shorten the working and setting times of the cement. Under warm conditions or if a long working time is required, it is therefore recommended that the glass slab be cooled. Testing to ISO 9917:1991 is carried out at 23 °C ± 1 °C and a relative humidity of 50 ± 10%.

³ Net setting time determined according to ISO 9917:1991.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

Manufacturer:**DENTSPLY DeTrey GmbH****De-Trey-Str. 1****78467 Konstanz****GERMANY****Phone +49 (0) 75 31 5 83-0****Distributor:****DENTSPLY Limited****Hamm Moor Lane****Addlestone, Weybridge****Surrey KT15 2SE****Phone (0 19 32) 85 34 22**

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© DENTSPLY DeTrey