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- b. Using the disposable brush supplied, immediately apply generous amounts of Prime & Bond® NT™ adhesive to thoroughly wet all the tooth surfaces. This surface should remain fully wet for 20 seconds and may necessitate additional applications of Prime & Bond® NT™ adhesive.
- c. Remove excess solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Surface should have a uniform glossy appearance. If not, repeat application and air dry.

- d. Cure adhesive for 10 seconds¹ using a curing light (i.e., ProLite™, Spectrum™).

- e. Apply a single coat of adhesive to the internal bonding surface of the restoration.Immediately air dry for 5 seconds.

Cementation Technique

- a. Dispense the desired shade of Calibra™ Esthetic Resin Cement base paste from the syringe directly onto the veneer. Protect cement from exposure to light.

- b. Seal the loaded veneer in place. Remove excess from the gingival margin with a blunt instru-ment.Brilliefly (10 seconds¹ or less) light cure the gingival priion only to tack restoration in place. Remove any excess from proximal and lingual margins.

Technique Tip: Placement of mylar strips between preparation and adjacent teeth prior to seating veneer aids in isolation and in excess cement clean up. After removing gingival excess and the systems in place, remove any remaining excess cement by pulling mylar strip towards the facial, i.e., from tooth to preparation.

- c. Light Cure 20 seconds¹ each from the buccal, lingual, and interproximal aspects.

Following the light cure check and adjust occlusion as necessary. Proceed to Finishing and Polishing, section 4.0.

2.0 CROWN/BRIDGE/INLAY/ONLAY CEMENTATION-DUAL CURED

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment

Application of Prime & Bond® NT™ Dual Cure Dental Adhesive System

- a. Place 1-2 drops of Prime & Bond® NT™ adhesive into a clean plastic mixing well. Replace cap promptly. Place an equal amount of drops of Self-Cure Activator into the same mixing well. Replace cap promptly. Mix contents for 1-2 seconds with a clean, unused brush tip.

- b. Using the disposable brush supplied, immediately apply mixed adhesive/activator to thoroughly wet all the tooth surfaces. These surfaces should remain fully wet for 20 seconds and may require additional applications of adhesive/activator.
- c. Remove excess solvent by gently drying with clean, dry air from a dental syringe for at least 5 seconds. Surface should have a uniform glossy appearance. If not, repeat application and air dry.

- d. Cure mixed adhesive/activator for 10 seconds¹ using a curing light (i.e., ProLite™, Spectrum™).

- e. Apply a single coat of mixed adhesive/activator to the internal bonding surface of the restoration. Immediately air dry for 5 seconds.

Cementation Technique

- a. Dispense the desired shade of Calibra™ Esthetic Resin Cement base paste from the syringe onto a clean mixing pad.Dispense an equal amount of the desired viscosity of catalyst paste (regular or High Viscosity Mix) into a clean mixing well.

Technique Tip: As with any dual cure resin cement system, the use of a dual cure adhesive system such as Prime & Bond® NT™ Dual Cure Dental Adhesive can shorten working time. The use of a dual cure adhesive system may be mixed with regular viscosity catalyst yielding a High Viscosity, “stay-put” rheotopic dual cured cement suitable for inlays, shallow inlays and veneers. Mixing Calibra™ Esthetic Resin Cement base with Regular Viscosity catalyst yields an easier flowing dual cured cement suitable for large, multi-surface and full coverage restorations with higher hydraulic seating pressure.

- b. Apply a uniform layer of cement on the entire internal surface of the restoration. For porcelain, it may be helpful to apply a thin layer of cement to the internal portions of the tooth preparation to avoid any porosity or voids. At room temperature, Calibra™ Esthetic Resin Cement, in either viscosity, works a work time of at least 2 min 30 sec.

- c. Seat the restoration with gradual pressure. A gentle rocking or vibratory motion may be helpful to insure optimal seating.

Remove gross excess from marginal areas. Use an instrument such as a blunted explorer, peri-odontal gross excess cleaner, dry brush tip. Restoration should not be moved or torqued during removal of gross excess cement.Special attention should be paid to interproximal areas using floss to remove excess cement.

Technique Tip: A 10 seconds¹ light “pre-cure” of excess cement at the margins will cause cement to “gel”, allowing easy cleanup.

Technique Tip: Apply moderate and consistent pressure to the restoration throughout the self-cure set time of approximately 6 minutes from the beginning of mixing.

- e. Once stabilized, light cure all marginal areas of the restoration for 20 seconds¹ from each direction - buccal, lingual and inter-occlusal.

Following the self-cure set check and adjust occlusion as necessary. Proceed to Finishing and Polishing, section 4.0.

3.0 Endodontic Posts

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment

Application of Prime & Bond® NT™ Dual Cure Dental Adhesive System

- a. Place 1-2 drops of Prime & Bond® NT™ adhesive into a clean plastic mixing well. Replace cap promptly. Place an equal number of drops of Self-Cure Activator into the same mixing well. Replace cap promptly. Mix contents for 1-2 seconds with a clean, unused brush tip.

- b. Apply mixed adhesive/activator to post preparation with the brush provided, being sure to apply generous amounts to the preparation orifices. A paper point pre-wetted with the adhesive mixture may be used in bringing the adhesive mixture down to the deepest portion of the preparation.Maintain contact of adhesive/activator with tooth structure for at least 20 seconds.

- c. Air dry treated post preparation with air syringe. Use of clean, dry paper points may be applied through removal of acetone in post space. Surface should have a uniform glossy appearance. If not, repeat application and air dry.

Following the self-cure set check and adjust occlusion as necessary. Proceed to Finishing and Polishing, section 4.0.

Cementation Technique

- a. Dispense the desired shade of Calibra™ Esthetic Resin Cement base paste from the syringe onto a clean mixing pad. Place an equal amount of Regular Viscosity catalyst paste. Mix the cement for 20-30 seconds. Dispense an equal amount of Regular Viscosity catalyst paste. Mix the cement for 20-30 seconds.

Technique Tip: As with any dual cure resin cement system, the use of a dual cure adhesive system such as Prime & Bond® NT™ Dual Cure Dental Adhesive can shorten working time.

- b. Spread mixed Calibra™ Esthetic Resin Cement components on surface of post and/or into the post preparation with a syringe tip, Lentulo Spinal or metal file.

- c. Seat post immediately. Clean up excess with appropriate instruments. A 10 seconds¹ light exposure “pre-cure” of excess cement at the margins will cause cement to “gel”, allowing easy cleanup.

- d. Stabilize post until cement sets. Apply moderate and consistent pressure to the restoration throughout the self-cure set time of approximately 6 minutes from the beginning of mixing.Light cure the cemented post for 20 seconds¹.

Proceed with core build-up and/or preparation.

4.0 Finishing and Polishing

Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of trauma to the enamel.

Polish final restoration using Prisma®-Gloss™ Polishing Paste and Prisma®-Gloss™ Extra Fine Polishing Paste with an Enhance® foam cup. Apply each paste with the foam cup first dry for 20-30 seconds then wet with increasing amounts of water in a wet polishing technique.

For technique assistance with the above-described DENTSPLY products call: 1-800-LD CALLK or 449-7531-583-333 (DENTSPLY DeTrey).

U.S. Patent Nos. 4,514,342;4,920,082;5,304,586 and other patents pending.

¹In Europe/DeTrey Conditioner 36

² Check curing light for minimum curing output of at least 300 mw/cm

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Die Entfernung von Zahnmaterialien ist am besten mit den Polierern, Ketten und Scheibchen des Enhance© Finerystystems zu erreichen. Das Enhance System entfernt die Überschüsse und poliert die Ränder der Restauration, ohne den Schmerz zu entfernen oder zu schädigen.