DENTSPLY
DETRAY

Dyract[®] Cem^{plus}

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Självhäftande compomercement

Self-Adhesive Compomer cementerings cement

Chemicky tuhnoúí fixaãni kompomer

Directions for Use _____	2
Gebrauchsanweisung _____	6
Mode d'emploi _____	10
Gebruiksaanwijzing _____	14
Istruzioni per l'uso _____	18
Instrucciones de uso _____	22
Instruções de utilização _____	26
Bruksanvisning _____	30
Brugsanvisning _____	34
Návod k použití _____	38

English

Deutsch

Français

Nederlands

Italiano

Español

Português

Svenska

Dansk

Česky

Dyract[®] Cem^{plus}

Self-Adhesive Compomer Luting Cement

Dyract[®] Cem^{plus} – a luting cement for dental use – combines the major benefits of glass ionomer cements – adhesion to tooth substance and fluoride release – with the mechanical strength of a luting composite.

Dyract Cem^{plus} is provided as a powder-liquid system in a translucent and an opaque version.

COMPOSITION

Powder:	Radiopaque strontium-alumino-silicate glass Titanium dioxide (opaque version only) Polymerisation initiator
Liquid:	Phosphate modified polymerisable monomers Carboxylic acid modified macromonomers Reactive diluent Polymerisation initiator Stabiliser

INDICATIONS

Cementation of conventional metal or porcelain-fused-to-metal inlays, onlays, crowns, bridges, posts and post-core units.

Adhesive cementation of porcelain or composite inlays, onlays and crowns¹. Cementation of Maryland bridges.

CONTRAINDICATIONS

1. Temporary luting of restorations.
2. Direct or indirect pulp capping:
Surfaces less than 1 mm from the pulp should have been covered with a hard-setting calcium hydroxide liner and be built up with a suitable material prior to final crown/inlay preparation.
3. Use on patients with a known allergy to methacrylates or other components of the cement.

WARNINGS

Dyract Cem^{plus} contains methacrylates, which may be irritating to skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of soap and water.

ADVERSE REACTIONS

The product may cause sensitisation by skin contact in susceptible persons. If skin sensitisation occurs discontinue use.

INTERACTIONS WITH DENTAL MATERIALS

Preferably use a non-eugenol containing temporary cement, as residual eugenol on the prepared tooth surface may interfere with setting and may cause softening of the polymeric components of the material.

¹ Instructions given by the manufacturer of the composite/ceramic material regarding indications, minimum crown thickness, materials and technique to be used for luting are to be observed. time of 20 seconds should be observed.

If a eugenol cement has been used, proper cleaning (as described below) is essential. If 3% H₂O₂ has been used, proper rinsing is essential. Higher concentration of H₂O₂ may interfere with the setting of the material and should not be used prior to the application of **Dyract Cem^{plus}**.

STEP-BY-STEP INSTRUCTIONS

1. Tooth and Restoration Pre-Treatment

1.1 Restoration

The fitting surfaces should be clean and dry prior to cementation.

Porcelain and composite: According to the manufacturers' recommendations, the fitting surface may be mechanically roughened (e.g. by sandblasting). In the case of porcelain, etching and silanisation may be required. Avoid touching the fitting surface afterwards.

Metal restorations: Sandblasting (50 µm aluminium oxide) of the fitting surface is recommended. In those instances where increased retention to noble metal is desired, additional tinplating may be used.

With **Dyract Cem^{plus}**, additional treatment of the restoration with an adhesive is not required.

1.2 Tooth

Clean prepared tooth surface and margins with non-oil containing, non-fluoridated pumice or prophylaxis paste (Nupro[®]) and a rubber cup on a low-speed handpiece. Wash thoroughly, isolate, and keep dry and uncontaminated by saliva, cervical fluid or blood.

2. Conditioning and Adhesive Application

In the case of standard retention requirements (comparable to glass ionomer cements), proceed to "3. Dosage and Mixing".

In the case of maximum retention requirements proceed to "2.1 Treatment of Preparation".

2.1 Treatment of Preparation

2.1.1 Conditioning with NRCTM 2

1. Dispense NRC into a DENTSPLY Applicator Dish³ or standard dappen dish.
2. Apply sufficient amounts of NRC with an Applicator Tip³ or disposable brush to prepared enamel and dentine. Leave undisturbed for 20 seconds. **Do not rinse.**
3. Remove excess NRC by blowing gently with an air syringe or blot dry with a cotton pellet. Do not desiccate the dentine structure.

2.1.2 Bonding with Prime&Bond[®] NT

One layer of Prime&Bond NT is applied:

1. Dispense Prime&Bond NT directly onto a fresh Applicator Tip³ or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish³ or standard dappen dish.
2. Immediately apply ample amounts of Prime&Bond NT to thoroughly wet all prepared tooth surfaces. This surface should be saturated which may necessitate additional application of Prime&Bond NT.
3. Leave the surface undisturbed for 20 seconds.
4. Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds. Surface should have a uniform, glossy appearance. If not, repeat steps 2 to 4.

² Alternatively, conventional acid-conditioning with a 34-36% phosphoric acid gel may be carried out: Gently extrude DeTrey[®] Conditioner 36 gel (36% phosphoric acid) to the cavity surfaces starting at the enamel margins. For best results, condition enamel for at least 15 seconds and dentine for 15 seconds or less. Remove gel with aspirator tube and/or vigorous water spray and rinse conditioned areas thoroughly for at least 15 seconds. Remove water from the rinsed cavity with a soft blow of air. Avoid desiccating the dentine, leave a moist surface.

³ DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

5. Light-cure for a minimum of 10 seconds⁴. Ensure uniform exposure of all prepared surfaces.

3. Dosage and Mixing

Shake the powder bottle to fluff the powder before dispensing. Overfill the scoop with powder and withdraw against the plastic leveller. Dispense appropriate quantities of powder and liquid onto a mixing pad. For dispensing the liquid hold bottle vertically over the mixing pad and gently squeeze the bottle. If the first drop contains air, disregard this drop.

The recommended mixing ratio is 1 scoop of powder to 1 drop of liquid. More powder (up to 2:1) may be added to achieve a more viscous mix.

2 scoops of powder and 2 drops of liquid will provide enough material for the luting of a well-fitting crown.

	Minutes
Mixing time (23 °C)	1/2
Extraoral manipulation time from start of mix (23 °C)	1 1/2
Intraoral working time (37 °C)	1
Setting time from insertion into oral cavity (37 °C)	4

Mix the material with a spatula until a homogeneous smooth paste is obtained. Mix within 30 seconds. Longer mixing and extraoral manipulation will shorten intraoral working and setting time. Shorter extraoral manipulation will prolong intraoral working and setting time.

Remove cement from mixing spatula and other instruments with a moist paper tissue before set of the material.

4. Luting

Line the restoration with **Dyract Cem^{plus}** and seat on the clean and dry tooth preparation. Do not apply cement into cavity or onto preparation.

Apply pressure and then release to allow the hydraulic pressure to dissipate. Repeat this procedure at least one more time and ensure that the restoration is completely seated.

Once in the mouth, the remaining working time is approximately 1 minute.

The material is completely set after approximately 4 minutes from insertion into the oral cavity.

5. Control of Excess Cement

Start immediately, while the material is still not set, to remove all excess cement.

This is best done with a plastic instrument and foam pellets applied with tweezers.

Interproximal removal is facilitated by the use of dental floss or Super Floss™ whilst holding the restoration stable with occlusal pressure.

Adhesion of excess cement to surrounding surfaces can be avoided by isolating teeth and restorations with Airblock™ gel prior to the application of cement.

Remove as much excess as possible before the cement sets.

Apply moderate pressure to stabilise the restoration until set of the cement has occurred (4 minutes).

If any excess cement remains after curing, remove with curette, scalers or rotary finishing burs.

⁴ When using a high-performance unit such as the curing lights manufactured by DENTSPLY, a curing time of 10 seconds is sufficient. For curing lights with an output lower than 300 mW/cm², a curing time of 20 seconds should be observed.

STORAGE

Keep out of sunlight.

To be stored at room temperature (23 °C). Storage at higher temperatures may reduce shelf stability. Storage of the material under refrigerated conditions will increase the viscosity of the liquid.

Effect of humidity: The Dyract Cem^{plus} powder and liquid bottles should be kept in a dry area. The powder bottle must be tightly closed immediately after use.

If the bottle has not been tightly closed, humidity may affect the setting of the material. In the case of water uptake, the desiccation cap will change its colour from blue to red.

BATCH NUMBER AND EXPIRY DATE

The batch number should be quoted in all correspondence that requires identification of the product.

Do not use after expiry date.

SPECIAL NOTES

Patent WO 96-19179.

Dyract Cem^{plus} contains an efficient initiator system which accelerates at mouth temperature (snap set reaction). Therefore, consider the following additional notes:

1. Use of Dyract Cem^{plus} at Elevated Temperatures

In case the temperature is higher than 23 °C, the setting reaction of **Dyract Cem^{plus}** is accelerated. Please consider the reduced working time.

2. Set of Material on the Mixing Pad

The hardening of **Dyract Cem^{plus}** is strongly temperature-dependant. Residual material on the mixing pad will at room temperature (23 °C) take approximately 6 minutes or longer to set and will exhibit a thick unpolymerised surface layer. Therefore, it cannot be used to judge the setting of the material in the mouth.

3. Air-Inhibited Layer

Under the conditions of the oral cavity, the air-inhibited layer on **Dyract Cem^{plus}** is less than 10 microns thick. For prevention of this layer on the surface of the cementation, Airblock gel may be applied on the crown/inlay margins after excess cement has been removed. Airblock application is, however, not mandatory.

For NRC and Prime&Bond NT please consult the respective Directions for Use for these products.

If you have any questions, please contact:

Manufacturer:
DENTSPLY DeTrey GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Hamm Moor Lane
Addlestone, Weybridge
Surrey KT15 2SE
Phone (0 19 32) 85 34 22

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© DENTSPLY DeTrey 1999-05-26