**DENTSPLY**  
DETREY

*Dyract<sup>®</sup>Seal*

**Compomer Pit & Fissure Sealant**

**Kompomer-Fissurenversiegler**

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***Dyract<sup>®</sup>Seal***

**Compomer Pit & Fissure Sealant**

Dyract® Seal is a light-curing dental pit & fissure sealant containing releasable fluoride. The fluoride offers an additional protection of the sealed pits and fissures.

Dyract Seal - when used in combination with NRC™ Non-Rinse Conditioner – offers a quick and simple application technique, which is of particular advantage with treatment of small children.

The application of Prime&Bond® NT prior to Dyract Seal allows a complete obliteration even of deep and narrow fissures (Total Seal™ Technique).

### **COMPOSITION**

Strontium-alumino-fluoro-silicate glass

Ammonium salt of PENTA (dipentaerythritol penta acrylate monophosphate) and N,N-dimethyl aminoethyl methacrylate

Carboxylic acid modified macromonomers

Diethyleneglycol dimethacrylate

Polymerisation initiator

Butylhydroxytoluene

Titanium dioxide (white version only)

### **INDICATIONS**

Preventive sealing of pits and fissures in the primary and secondary dentition.

### **CONTRAINDICATIONS**

Known allergy to dimethacrylate resins and other components of the material.

Not to be used in case of caries.

### **WARNINGS/PRECAUTIONS**

Avoid contact of Dyract Seal with oral tissues, eyes and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with soap and water.

Dyract Seal may cause sensitisation by skin contact in susceptible persons. If sensitisation occurs, discontinue use. Remove soft (oxygen-inhibited) surface layer after light-curing with cotton pellets or cotton rolls.

## **ADVERSE REACTIONS**

Reversible inflammatory changes of the oral mucosa after accidental contact have been associated with (meth)acrylate monomers.

## **STEP-BY-STEP INSTRUCTIONS**

### **Cleaning**

Clean surfaces to be sealed with oil free paste such as Nupro® Prophy Paste. Rinse well with water.

### **Isolating**

Isolate the teeth to be sealed with rubber dam or cotton rolls and dry each tooth with air, free of oil or water contamination.

Dyract Seal may be used in two different application techniques:

## **A. Conventional Technique**

### **1. Application of Phosphoric Acid Conditioner**

Apply 35 - 37% phosphoric acid liquid or gel (e. g. DeTrey® Conditioner 36) to the fissures and condition enamel for at least 30 seconds.

### **2. Rinsing and Drying**

Remove acid liquid/gel with aspirator tube and/or vigorous water spray and rinse for at least 15 seconds. Dry surface with a blow of air. Properly etched enamel has a dull, frosty white appearance. If this is not the case, re-etch enamel, rinse and dry as described above.

Once the fissure has been properly treated, it must be kept uncontaminated. If saliva contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure for 5 seconds only. Rinse and dry as described above.

## **B. Total Seal Technique**

### **1. Application of NRC Non-Rinse Conditioner**

- 1.1 Dispense NRC into a DENTSPLY Applicator Dish<sup>1</sup> or standard dappen dish.
- 1.2 Apply 1 drop of NRC with an Applicator Tip<sup>1</sup> or disposable brush to the fissure. Leave undisturbed for 20 seconds. **Do not rinse.**
- 1.3 Remove excess NRC by blowing gently with an air syringe.

Once the fissure has been properly treated, it must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with forceful water-spray and repeat the application of NRC.

### **2. Application of Prime&Bond NT**

- 2.1 Dispense Prime&Bond NT directly onto a fresh Applicator Tip<sup>1</sup> or onto a disposable brush. Alternatively, dispense into a fresh DENTSPLY Applicator Dish<sup>1</sup> or standard dappen dish.

- 2.2 Immediately apply 1 drop of Prime&Bond NT to the fissure.
- 2.3 Leave undisturbed for 20 seconds.
- 2.4 Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds.

### **Application of Dyract Seal**

1. Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto pad.
2. Attach disposable needle tip to end of the syringe. **Turn tip clockwise 1/4 to 1/2 turn to assure that it is fully seated.** Tug on tip to be sure that it is locked into the collar of the syringe.
3. Dyract Seal should flow freely with gentle pressure. Do not use excessive force. If more than gentle pressure is required, remove from patient field and check for obstruction.
4. Dispense Dyract Seal directly into the fissure.
5. Discard needle immediately after use. Replace original cap. Do not store syringe with dispensing tip in place. Store only with original cap.

---

<sup>1</sup> DENTSPLY Applicator Dish and Applicator Tips are available from your dental dealer.

6. It is recommended to pull back slightly on the syringe plunger after use to prevent excessive flow of material.

### **Curing**

Cure for at least 20 seconds with a dental polymerisation unit (e. g. DENTSPLY® QHL 75™, Spectrum® 800) keeping the tip of the light guide as near as possible to the tooth without touching it.

Remove soft (oxygen-inhibited) surface layer after light-curing with cotton pellets or cotton rolls.

### **Retention and Occlusion Control**

Check with an explorer for complete coverage and retention. Check with articulating paper for premature occlusal contacts and correct, if necessary, with a finishing bur.

### **STORAGE**

The Dyract Seal syringes should be tightly closed immediately after use.

Keep out of direct sunlight.

Do not store at temperatures exceeding 24 °C.

Humidity can adversely affect the properties of unsealed Dyract Seal syringes. Therefore keep syringes sealed in their blister pack until first use. Under high humidity (90% r.h.) unsealed syringes stay usable for at least 3 weeks.

Under normal ambient conditions, unsealed syringes stay usable for about 12 weeks. Discard syringes that require more than gentle pressure to extrude Dyract Seal.

### **BATCH NUMBER AND EXPIRY DATE**

The batch number should be quoted in all correspondence which requires identification of the product.

Do not use after expiry date.

If you have any questions, please contact:

**Manufacturer:**

**DENTSPLY DeTrey GmbH**  
**De-Trey-Str. 1**  
**78467 Konstanz**  
**GERMANY**  
**Phone +49 (0) 75 31 5 83-0**

**Distributor:**

**DENTSPLY Limited**  
**Hamm Moor Lane**  
**Addlestone, Weybridge**  
**Surrey KT15 2SE**  
**Phone (0 19 32) 85 34 22**

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© DENTSPLY DeTrey 2000-12-18