Directions for use

# Enhance®

Composite Finishing & Polishing System

## INDICATIONS

The Enhance® Composite Finishing & Polishing System is designed for use in the finishing of anterior or and posterior composites or any other area that requires finishing. When used as directed, this system will yield a highly polished surface on hybrid composites such as TPH® Spectrum™, Prisma® TPH®, and Prisma® APH® Composites and SureFill™ High Density Restorative Material as well as microfilled composite restoratives. The Enhance® Finishing & Polishing System is useful for finishing and polishing compomers such as Dyract® AP.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a dentist.

## COMPOSITION

The Enhance® Finishing & Polishing System is comprised of the Enhance® Finishing Disc, Finishing Cup, Finishing Point, the Enhance® Polishing Cup, the Enhance® Cup Mandrel, Prisma®-Gloss™ and Prisma®-Gloss™ Extrafine Polishing Pastes.

## CONTRAINDICATIONS

- Enhance® Finishing & Polishing System Discs, Cups and Points are recommended as a single use disposable item. Attempts to sterilize these items can adversely affect their performance.
- ENHANCE® FINISHING DISCS, POINTS AND CUPS SHOULD NOT BE USED IN HIGH SPEED HANDPIECES.

## PRECAUTIONS

- If polishing a microfilled composite restoration, polishing is complete after Prisma®-Gloss™ only. It is not necessary to use the Prisma®-Gloss™ Extrafine.
- As with any rotary instrument, heat will build up when using the Enhance® Finishers with prolonged contact. Use with intermittent pressure.
- Use flat end and corner edge of Enhance® Polishing Cup, excessive lateral pressure may dislodge cup from mandrel.
- Care should be taken to have good control of any rotary instrument to protect patients from injury.

## STEP BY STEP INSTRUCTION FOR USE

**STEP 1:** Place restoration or composite veneer as per standard techniques. Complete gross reduction of excess and general outline form of the restoration using Prisma® Finishing Burs, or diamond finishing instrument.

**STEP 2:** Insert an Enhance® Finisher (Disc, Point or Cup) into a conventional speed contra-angled handpiece and continue finishing. **NOTE: ENHANCE® DISCS, POINTS AND CUPS MUST NOT BE USED IN HIGH SPEED HANDPIECES.** The aggressiveness of the Enhance® Finishers is controlled by the pressure applied to the surface of the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.

**STEP 3:** Attach an Enhance® Polishing Cup to the provided mandrel by inserting the bayonet end of the mandrel into the opening on the narrow end of the cup. Rotate the mandrel 1/4 turn makes insertion easier. Make sure the mandrel inserts fully into the cup. Insert the mandrel into a conventional speed contra-angles handpiece.

**STEP 4:** Apply a small amount of Prisma®-Gloss™ Paste to the surface of the Enhance® Polishing Cup. Work the surface of the restoration, dry initially, at moderate speed and pressure.

**STEP 5:** To increase surface luster, ADD WATER IN SMALL AMOUNTS (i.e. dropwise) to dilute the paste using a light circular buffing motion. Repeat as needed to produce a smooth surface.

**STEP 6:** Rinse Prisma®-Gloss™ Paste from the tooth surface and Enhance® Polishing Cup. **THIS IS THE FINAL STEP FOR MICROFILLED COMPOSITES.** Continue with the next step for hybrid composite materials like TPH® Spectrum™, Prisma® TPH®, and Prisma®-APH® Composite and Dyract® AP Compomer.

**STEP 7:** Apply Prisma®-Gloss™ Extrafine Paste to the Enhance® Polishing Cup as described in step 4 and lightly buff the surface for the final luster. Add water and continue polishing. A high luster should rapidly develop in 20-30 seconds of polishing. Rinse all surfaces and discard the Enhance® Cup after use.

## TECHNIQUE HINTS:

- The aggressiveness of the Enhance® Disc is controlled by the pressure applied to the surface of the composite. The greater the pressure, the more material is removed; lighter pressure leaves a smooth surface without removing bulk.
- Enhance® Finishing Points are ideal for posterior occlusal and concave lingual surfaces. Enhance® Finishing Cups are well suited to cervical and proximal line angle areas.
- Polishing efficiency is greatly increased and splatter of Prisma®-Gloss™ is greatly reduced by rubbing paste into the surface of the Enhance® Polishing Cup before use. This procedure leaves a thin, adherent surface film of paste on the cup surface. Use the paste impregnated cup dry at low speed; after 30-60 seconds, apply a small amount of water to the tooth and cup. Use again at low speed for another 30-60 seconds to achieve high surface luster.

©1998 DENTSPLY International Inc. All Rights Reserved  
Manufactured by  
DENTSPLY Caulk  
DENTSPLY International Inc.  
Milford, DE 19963-0359  
U.S.A.

Authorized EU Representative  
DENTSPLY DeTrey GmbH  
D-78467 Konstanz  
Germany

Swiss Representative  
DENTSPLY DeTrey GmbH, Konstanz  
Zweigiederlassung Baar  
Oberdorfstrasse 11  
CH 6342 Baar

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