

x·smart[®]plus

everything you like about x·smart[®] with a plus



now with
PROTAPER·NEXT™
settings

+
**WE
KNOW
ENDO.**

everything you like about x•smart[®] with a plus



Pre-programmed for different file systems.
Choose your preferred system at a single glance



User friendly interface with
«click and go» navigation
and button controls dedicated
to each function

Excellent visibility and access
due to the miniature contra-angle

No foot pedal
On/Off button on the
motor handpiece



Reciprocation

+ Single file root canal shaping
with the wave•one[®] reciprocating motion



Continuous motion

+ The next generation of the endo gold standard :
PROTAPER•NEXT™, the performance of PROTAPER[®] UNIVERSAL
with added versatility and a shorter clinical sequence



+ The file library also includes preset programs for PROTAPER[®] UNIVERSAL, PATHFILE[®], GATES, RECIPROC, as well as 7 free programs for individual settings in continuous rotation

technical features

- Miniature contra-angle head that can be adjusted in 6 positions
- On/Off button on the motor handpiece
- Reciprocating motion and continuous rotation
- Includes a file library with preset programs for wave•one®, PROTAPER•NEXT™, PROTAPER® UNIVERSAL, PATHFILE®, GATES, RECIPROC – 7 free programs for individual settings in continuous rotation
- Auto reverse rotation at preset torque limit
- Torque range : 0.6 – 4.0 Ncm
- Speed range : 250 – 1'200 rpm
- Operates on a rechargeable battery and when connected to the mains
 - Battery autonomy : approx. 2 hours continuous use
 - Full recharge in approx. 5 hours



X-smart^{plus}

Endo motor with contra-angle and universal AC adaptor
REF A 1032



X-smart^{plus} 6:1 contra-angle

REF A 1033



X-smart^{plus} cleaning adaptor

Cleaning and maintenance unit adaptor
REF A 1007-600